



離心濃縮機附件

ACCESSORIES FOR CES-8000, CES-8080

ROTORS

ROTOR for CES-8000

R-8000A, 1.5ml x 12 pcs

R-8000B, 1.5ml x 24 pcs

R-8000C, $\phi 13 \times 100\text{mmL} \times 8$ pcs or $\phi 12 \times 57\text{mmL} \times 16$ pcs

Specified for scintillation vial is available



ROTOR for CES-8080

R-8080A, 1.5ml x 72 pcs

R-8080B, 96-well plate x 2 sets


R-8080C, 15ml falcon tube x 20 pcs

R-8080D, 50ml falcon tube x 8 pcs

R-8080E, $\phi 12 \times 75\text{mmL} \times 48$ pcs

R-8080F, $\phi 13 \times 100\text{mmL} \times 48$ pcs

Specified for other tube or scintillation vial is available.



COLDTRAP



MODEL: CT-5020, - 50°C

MODEL: CT-5000, - 50°C

MODEL: CT-8500, - 85°C

CT-5020

VACUUM PUMP



G-1000

GVD-050A, 60ℓ/min.

G-100D, 120ℓ/min.

GLD-201B, 240ℓ/min.

Vacuum: $6.7 \times 10^{-2}\text{Pa}$

Optional accessory:
Oil Mist Trap

組合例



DNA CONCENTRATOR



CENTRIFUGAL EVAPORATOR
COMPLETE SET



CENTRIFUGAL EVAPORATOR &
FREEZE DRYER COMPLEX SET