



## 經濟型真空冷凍乾燥機

ECONOMIC FREEZE DRYER CT-5020D, CT-5000D, CT-8500D



CT-5020D+MF-4

CT-5000D

CT-8500D

\*The above Manifold (MF-4) is optional accessory  
Other usages, please refer to Freeze dryer accessories.

### SPECIFICATIONS

MODEL	CT-5020D	CT-5000D	CT-8500D
TRAP TEMP. *	-50°C		-85°C
TRAP CAPACITY	3.5ℓ (1ℓ/batch)	4.5ℓ (1.5ℓ/batch)	4.5ℓ (1.5ℓ/batch)
TRAP MATERIAL	SUS304 plus Teflon coating		
TEMP. DISPLAY	Optional		
COMPRESSOR	1		2
REFRIGERANT	CFC-free		
DEFROST	Hot water		
SAFETY DEVICE	Power leakage, over current prevention devices		
WEIGHT	36 kg	40 kg	64 kg
DIMENSIONS (WxDxH)	340x450x490 mm	320x420x685 mm	360x500x840 mm
POWER SOURCE	110/220V, 50/60Hz		220V, 50/60Hz
OPTIONAL ACCESSORIES	Manifold, Drying chamber (or plus sealing device) Vacuum pump, Drying flask		

※At R.T.20°C

\*\* (本公司保有修改規格之權利, 不另行通知)

We reserve the rights to change the specifications without notice.)

### 特點:

- 經濟的價格。
- 可加裝LED溫度顯示器。
- 使用熱水除霜方式, 並附"取冰手", 可輕易將冰取出。
- 方便使用。