



真空冷凍乾燥機

FREEZE DRYER, BENCH TOP FD-5020



FD-5020+MF-4

FD-5020+MF-8

FD-5020+MF-18

FD-5020+DC-1000

*The above Manifold (MF-8, MF-18, MF-4) & Drying chamber (DC-1000) are optional accessories.
Other usages, please refer to Freeze dryer accessories.

SPECIFICATIONS

MODEL	FD-5020
TRAP TEMP. °	-54°C
TRAP CAPACITY	3.5L (1L/batch)
TRAP MATERIAL	SUS304 plus Teflon coating
TEMP. DISPLAY	LED display
TEMP. DISPLAY RANGE	-100 ~ +100°C
VACUUM DISPLAY	LED display
VACUUM DISPLAY RANGE	1x10 ³ ~ 1x10 ⁻¹ Pa
VACUUM GAUGE	Pirani vacuum gauge
COMPRESSOR	1
REFRIGERANT	CFC-free
DEFROST	Heater
SAFETY DEVICE	Power leakage, over current prevention devices
WEIGHT	38 kg
DIMENSIONS (WxDxH)	340 X 450 x 490 mm
POWER SOURCE	110/220V, 50/60Hz
OPTIONAL ACCESSORIES	Manifold, Drying chamber(or plus sealing device) Vacuum pump, Drying flask

※At R.T.20°C

** (本公司保有修改規格之權利, 不另行通知。)

We reserve the rights to change the specifications without notice.)