



# 真空冷凍乾燥機

FREEZE DRYER FD-5030, FD-5060, FD-8530, FD-8560



FD-5030+MF-8

FD-8530+MF-18

FD-8560/FD-5060+DC-1000+MF-4

\*The above Manifold (MF-8, MF-18, MF-4) & Drying chamber (DC-1000) are optional accessories.

## SPECIFICATIONS

MODEL	FD-5030	FD-5060	FD-8530	FD-8560
TRAP TEMP. °	-54°C		-85°C*	
TRAP CAPACITY	4.5ℓ (1.5ℓ/batch)	8ℓ (3ℓ/batch)	4.5ℓ (1.5ℓ/batch)	8ℓ (3ℓ/batch)
TRAP MATERIAL	SUS304 plus Teflon coating			
TEMP. DISPLAY	LED display			
TEMP. DISPLAY RANGE	-100 ~ +100°C			
VACUUM DISPLAY	LED display			
VACUUM DISPLAY RANGE	1x10 <sup>-3</sup> ~ 1x10 <sup>-1</sup> Pa			
VACUUM GAUGE	Pirani vacuum gauge			
COMPRESSOR	1		2	
REFRIGERANT	CFC-free			
DEFROST	Heater			
SAFETY DEVICE	Power leakage, over current prevention devices			
WEIGHT	42 kg	50 kg	66 kg	68 kg
DIMENSIONS (WxDxH)	320x420x685 mm		360 x 500 x 840 mm	
POWER SOURCE	110/220V, 50/60Hz		220V, 50/60Hz	
OPTIONAL ACCESSORIES	Manifold, Drying chamber(or plus sealing device) Vacuum pump, Drying flask			

※At R.T.20°C

\*\* (本公司保有修改規格之權利, 不另行通知)

We reserve the rights to change the specifications without notice.)