



## 真空冷凍乾燥機附件

### ACCESSORIES FOR FREEZE DRYER

#### 多歧管 MANIFOLD

#### SPECIFICATIONS



COLUMN	φ 76 x 360mm L	φ 76 x 320mm L	φ 76 x 210mm L
VALVE Number	18 pcs	8 pcs	4 pcs
VALVE	Rubber, 10mm <sup>1</sup>		
BASE	Transparent		
OPTIONAL ACCESSORY	Multi-connector		

#### 圓筒乾燥室 DRYING CHAMBER



DC-1000

1. DC-1000
  - Single plate
  - φ250 x 250mm L
  - 30°C ± 2°C
2. DC-1000S  
DC-1000 plus sealing device
3. DC-3000
  - Triple plate
  - φ250 x 380mm L
  - 30°C ± 2°C
4. DC-3000S  
DC-3000 plus sealing device

※DC系列可選配加裝溫度控制 (溫度範圍: 30~42°C)  
DE系列為無加熱式

#### 廣口軟蓋乾燥瓶 DRYING FLASK



- 150ml (φ 60mm)
- 300ml (φ 60mm)
- 600ml (φ 90mm)
- 1000ml (φ 90mm)
- 1200ml (φ 90mm)
- 2000ml (φ 90mm)

平底, 瓶身附刻度, 矽膠製軟蓋, 標準聯接管為 3/4".  
※可更換聯接管, 能適合任一品牌的冷凍乾燥機