



## 轉動預凍槽

PREFREEZING BATH PFB-5000



PFB-5000

### 特點

- 使液體凍結於瓶壁, 冷凍乾燥過程時間縮短並避免凍乾瓶於冷凍櫃中預凍時因樣品過量導致破裂.
- 適用最大1000ml標準凍乾瓶, 及磨砂 29或 24 2000ml的梨形瓶.
- 馬達角度調整: 可依樣品容量不同, 快速拉鈕調整三種角度
- 轉動軸心長度可調
- 內藏抽氣幫浦可抽轉動瓶
- 具自動洩壓裝置與抽氣幫浦連動

凍乾瓶連接樣式



梨形型瓶連接樣式



### SPECIFICATIONS

MODEL	PFB-5000
BATH TEMP. *	-54°C
BATH CAPACITY	6L
BATH MATERIAL	SUS304
TEMP. DISPLAY	LED display
TEMP. DISPLAY RANGE	-100 ~ +100°C
VACUUM PUMP	Build-In
RELEASE VALVE	Build-In, Electric
SPEED CONTROL	10~160rpm
COMPRESSOR	1
REFRIGERANT	CFC-free
SAFETY DEVICE	Power leakage, over current prevention devices
DIMENSIONS (WxDxH)	360 x 500 x 1000 mm
POWER SOURCE	110/220V, 50/60Hz

※At R.T. 20°C

\*\* (本公司保有修改規格之權利, 不另行通知)

We reserve the rights to change the specifications without notice.)