



# 小型真空烘箱(方型)

## MICRO VACUUM OVEN VO-2000



VO-2000

### 用途

- 化學實驗之溶媒及水份抽乾.
- 分析實驗, 樣品之製備.
- 生化實驗, 生化樣品微量水份抽乾使用.
- 材料實驗室, 樣品乾燥使用.
- 各式微量樣品乾燥使用.

### 特點

- 少量樣品乾燥最適宜, 適用於各式閃爍計數瓶, 試管, 培養皿, 小燒瓶, 小燒杯等樣品容器.
- 透明式上蓋, 易於觀察乾燥過程.
- 不需要繁雜的操作步驟.
- 體積小, 不暫實驗室空間.
- 直結式抽氣, 洩壓接頭, 無銅管及軟管連接, 降低污染.

### SPECIFICATIONS

MODEL	VO-2000
TEMP. RANGE	Room temp.~200°C
TEMP. ACCURACY	±1°C
INNER DIMENSIONS(WxDxH)	200 x 300 x 200mm
TOTAL CAPACITY	12ℓ
TEMP. CONTROL	PID control
ULTIMATE VACUUM	< 1 Torr
VACUUM GAUGE	0~760mmHg
HEATING TYPE	Direct heating
HEATER	600W
VISIBALE WINDOW	8mm, tempered glass
SHLF	1pc, stainless steel
TOTAL DIMENSIONS(WxDxH)	300 x 380 x 420mm
POWER SOURCE	110/220V, 50/60Hz, 5A

(本公司保有修改規格之權利, 不另行通知.)

We reserve the rights to change the specifications without notice.)

### 組配圖及其選擇性配件:



1. 雙層組合架(RACK 60)
2. 低溫冷凝槽 Cold Trap (免加抗凍劑型)  
CT-5000 -50°C  
CT-8500 -85°C
3. 真空幫浦 GVD-050A  
60ℓ/min., 5x10<sup>4</sup> Torr