



離心濃縮機

CENTRIFUGAL EVAPORATOR CES-8000, CES-8080



特點:

- Gentle concentration and drying under non-bumping and non-foaming.
- 可選用配件: 真空幫浦, Cold Trap 及 轉子.

SPECIFICATIONS

MODEL	CES-8000	CES-8080
ROTOR SPEED	100~1,800 rpm with meter	
TEMP. RANGE	R. T. +5°C~80°C	
TEMP. STABILITY	±0.1°C	
TEMP. SETTING & DISPLAY	LED display	
DRIVE SYSTEM	Magnetically coupled	
ULTIMATE VACUUM	Bellow 10 ⁻¹ Torr	
VACUUM GAUGE	0~760mmHg	
CHAMBER MATERIAL	SUS304 plus Teflon coating	
CHAMBER DIMENSIONS	φ142 x 110mm H	φ255 x 165mm H
ROTOR DIMENSIONS	Optional	
SAFETY DEVICES	Over current prevention devices	
WEIGHT	14 kg	25 kg
TOTAL DIMENSIONS (WxDxH)	330 x 300 x 290mm	400 x 350 x 380mm
POWER SOURCE	110/220V, 50/60Hz	

** (本公司保有修改規格之權利, 不另行通知。)

We reserve the rights to change the specifications without notice.)