



乾式真空幫浦, 耐腐蝕型

DIAPHRAGM VACUUM PUMP N838.3, N820.3, N840.3, N920 / KNF



用途:

- 特別適用於旋轉真空濃縮機.
- N838.3及N920適用於中度酸鹼溶媒環境的濃縮乾燥.
- N820.3及N840.3適用於高度腐蝕性的氣體及液體.
- 全系列都採用TEFLON製的膜片.
- 出口側可另選配廢液回收裝置.

SPECIFICATIONS

MODEL	N838.3	N820.3	N840.3	N920
PUMPING SPEED	22ℓ/min.	20ℓ/min.	34ℓ/min.	20ℓ/min.
ULTIMATE PRESSURE	15mbar	8mbar	8mbar	2mbar
PUMP HEAD	PPS	PTFE	PTFE	PPS
DIAPHRAGM	PTFE-coated			
VALVE	FFPM			
INLET PORT DIAMETER	φ 10 mm	φ 10 mm	φ 10 mm	φ 9 mm
OPERATING CURRENT	0.5A	0.7A	1.5A	1.4A
MOTOR	80W	120W	245W	135W
MOTOR PROTECTION	IP20	IP44	IP44	IP20
TOTAL DIMENSIONS(WxDxH)	404x110x210 mm	312x154x207 mm	341x166x226 mm	324x158x226 mm
WEIGHT	6.8 kg	9.3 kg	12.6 kg	8.6 kg
POWER SOURCE	110V, 60Hz			

(本公司保有修改規格之權利, 不另行通知。
We reserve the rights to change the specifications without notice.)