



液態氮容器

LIQUID NITROGEN CONTAINER

◎運輸及儲存液態氮桶, 適合實驗室工作場所

SPECIFICATIONS

MODEL	LNC-10	LNC-15	LNC-30	LNC-35	LNC-50
CAPACITY (ℓ)	10	15	30	35	50
NECK DIA. (mm)	50				
OUTER DIA. (mm)	294	354	445	445	445
HEIGHT (mm)	525	600	660	690	840
WEIGHT EMPTY (kg)	6	8.5	13.5	15.5	17
WEIGHT FULL (kg)	14	20.6	37.7	43.8	57.4
HOLDING TIME (day)	≥42	≥65	≥127	≥130	≥165
PROTECTIVE LEATHER	Standard				
WHEELS	—				Standard

LNC-30/LNC-35/LNC-50選用配備:

- 液態氮泵 LNP-4, 具閥門, 壓力錶, 安全閥
- 全系列可選用Canister, 可做為生物樣品儲存, 變成兩用桶.

◎生物樣品保存, 超低蒸發率液態氮桶

SPECIFICATIONS

MODEL	LNC-3B	LNC-6B	LNC-10B	LNC-15B	LNC-30B	LNC-35B	LNC-80B	LNC-100B
CAPACITY (ℓ)	3	6	10	15	30	35	80	100
NECK DIA. (mm)	50					80	200	
OUTER DIA. (mm)	234	294	294	354	445	445	556	556
HEIGHT (mm)	400	420	525	595	660	690	980	1100
WEIGHT EMPTY (kg)	3	5	5.6	8	13	15	28	33
WEIGHT FULL (kg)	5.5	9.8	13.7	20.1	37.2	43.3	93	114
HOLDING TIME(day)	≥22	≥42	≥70	≥112	≥220	≥140	≥80	≥100
NO. OF CANISTER	3			6			3	
PROTECTIVE LEATHER	Standard							
WHEELS	—						Standard	

◎生物樣品保存, RACK型式

SPECIFICATIONS

MODEL	LNC-35R	LNC-100R
CAPACITY (ℓ)	35	100
NECK DIA. (mm)	125	290
OUTER DIA. (mm)	445	556
HEIGHT (mm)	690	1100
WEIGHT EMPTY (kg)	16	33
WEIGHT FULL (kg)	44.3	114
HOLDING TIME(day)	≥95	≥70
RACK TYPE CANISTER	6	4
RACK/CANISTER	5	8
2.0ml VIAL/RACK	25	81
TOTAL 2.0ml vial (pcs)	750	2592



(本公司保有修改規格之權利, 不另行通知。)

We reserve the rights to change the specifications without notice.)