



旋轉真空濃縮機 — 直立式, 1ℓ

ROTARY EVAPORATOR R-2000V, R-2000V-B1, R-2000V-B2



SPECIFICATIONS

MODEL	R-2000V	R-2000V-B1	R-2000V-B2
CONDENSER		Vertical	
SAMPLE FLASK		Standard 1ℓ (T29/32), Max. 3ℓ	
RECEIVER FLASK		1ℓ (S35/20)	
ROTATION SPEED		10~200 rpm	
ULTIMATE VACUUM		Bellow 1 Torr	
ROTARY JOINT (φ x L)		φ25 x 185 mm, T29/32	
JACK HEIGHT		Main unit, 110 mm	
TEMP. RANGE		30°C ~ 100°C	
TEMP. SET/DISPLAY	—	Analog	LED (0.5°C)
TEMP. ACCURACY	—		±1.5°C
BATH CAPACITY(φ x H)	—	φ275(φ210, bottom) x 110mm	
TOTAL DIMENSIONS(WxDxH)	520x270x820mm	570 x 340 x 820 mm	
WEIGHT	12.1 kg	14.2 kg	14.4 kg
POWER SOURCE		110V/220V, 50/60Hz	



濃縮瓶浸入角度調整
電子控制主體是正的
(R-2000/R-3000系列共通)

(本公司保有修改規格之權利, 不另行通知。
We reserve the rights to change the specifications without notice.)